



Email : sales@naugramedical.com  
Phone: 9896600003

**Product Name :**  
Miller Type Straight Laryngoscope Blades - Stainless Steel  
(Polished Finish)

**Product Code :**  
NMLYG10A08



**Description :**

Miller Type Straight Laryngoscope Blades - Stainless Steel (Polished Finish) Laryngoscope Blade SS. Blade Size 3. Polished. Autoclavable

**Technical Specification :**

The most popular of the straight blade types.

When originally designed, the size of the flange was reduced to minimize trauma and the curve at the beak or top was extended to improve lifting of the epiglottis.

These improvements facilitate greater exposure of the larynx in difficult patients.

This blade is available in a variety of sizes.

Blade sizes 2, 3, and 4 have lamp on right side. Blades made of surgical stainless steel.



**Medical Hospital Equipments Manufacturers**

Website: Medical Hospital Equipments Manufacturers,

