

Product Name :
Macintosh Type Disposable Fibre Optic Laryngoscope
Blade (Plastic heel, SS Spatula)

Product Code :
NMLYG10D02



Description :

Macintosh Type Disposable Fibre Optic Laryngoscope Blade (Plastic heel, SS Spatula) Laryngoscope Blade - Plastic heel, SS Spatula

Technical Specification :

Laryngoscope set is used for larynx examination, anesthesia, and resuscitation for endotracheal intubation.

Macintosh type curved blades are available in dull finish hence reduce reflections. Plastic handles use type C battery that lasts longer.

Laryngoscopes are simply specialized torches with a handle containing the batteries and a blade which holds the bulb.