



Email : sales@naugramedical.com  
Phone: 9896600003

**Product Name :**  
Miller Type Straight Laryngoscope Blades - Stainless Steel  
(Satin/Dull Finish)

**Product Code :**  
NMLYG10A17



**Description :**

Miller Type Straight Laryngoscope Blades - Stainless Steel (Satin/Dull Finish) Laryngoscope Blade SS. Blade Size 4. Matt

**Technical Specification :**

Combines the straight Jackson blade design with the curved distal tip Miller design to allow greater visibility and a 'direct line' approach to the trachea  
Low profile flange reduces risk of damage to teeth  
Crafted in surgical stainless steel  
Satin matte finishes which eliminates reflection and glare  
Fit any Standard Conventional handles



**Medical Hospital Equipments Manufacturers**

**Website:**Medical Hospital Equipments Manufacturers,