



Email : sales@naugramedical.com
Phone: 9896600003

Product Name :
Macintosh Type Reusable Fibre Optic Laryngoscope Blade

Product Code :
NMLYG10B16



Description :

Macintosh Type Reusable Fibre Optic Laryngoscope Blade Laryngoscope Blade - With in built Fibre Optic Tube Autoclavable

Technical Specification :

High grade of fibre-optics provides optimal light transfer without heat, reducing risk of thermal damage and maximising view of the larynx.

Made from high quality stainless steel, our blades give clinicians a complete choice for all patient intubation requirements.



Medical Hospital Equipments Manufacturers

Website: Medical Hospital Equipments Manufacturers,