

Product Name :
Battery Handle for Laryngoscope Blades

Product Code :
NMLYG10H05



Description :

Battery Handle for Laryngoscope Blades Size: Penlight. Type: Polished Finish. Material: Stainless Steel. Suitable for Laryngoscope Blades. Autoclavable

Technical Specification :

Laryngoscope handle offers many clear advantages, including a preferred cold-light technology, a halogen lamp inside the handle, a battery-powered handle made with high quality polycarbonate material.

It uses two 1.5 volt alkaline C-Cells or two 1.2 volt nickel metal-hydride rechargeable batteries.

NAUGRA[®]

Medical Hospital Equipments Manufacturers

Website: [MedicalHospitalEquipmentsManufacturers,](http://MedicalHospitalEquipmentsManufacturers.com)