



Email : [sales@naugramedical.com](mailto:sales@naugramedical.com)  
Phone: 9896600003

**Product Name :**  
Battery Handle for Laryngoscope Blades

**Product Code :**  
NMLYG10H06



**Description :**

Battery Handle for Laryngoscope Blades Size: Penlight. Type: Satin/Dull Finish. Material: Stainless Steel. Suitable for Laryngoscope Blades. Autoclavable

**Technical Specification :**

Laryngoscope handle offers many clear advantages, including a preferred cold-light technology, a halogen lamp inside the handle, a battery-powered handle made with high quality polycarbonate material.

It uses two 1.5 volt alkaline C-Cells or two 1.2 volt nickel metal-hydride rechargeable batteries.



**Medical Hospital Equipments Manufacturers**

Website: [SITE](http://SITE)