

**Product Name :**  
Laryngoscope Handle with Inbuilt Led Bulb**Product Code :**  
NMLYG10H07**Description :**

Laryngoscope Handle with Inbuilt Led Bulb Autoclavable

**Technical Specification :**

Light weight construction.  
Diamond knurled handle made surface a secure grip in the hand.  
Practically wear free, standard holding fixture for the blade.  
Solid handle made of chromium plated metal.  
Solid metal contact for reliable activation of illumination.  
Simple change of batteries at the base of the handle.