

Product Name :
Miller Laryngoscope Set (With Handle)**Product Code :**
NMLYG10S06**Description :**

Miller Laryngoscope Set (With Handle) The above set are satin finish all stainless steel. Other custom made combination of sets available upon request

Technical Specification :

The blades of the laryngoscope are made of 304-grade stainless steel. The material makes the blade safe to use and durable.

The handle is made of durable stainless steel. It sports tiny grooves all over its body to allow a good grip.

An LED lamp is fitted along the blade to help a trained medical practitioner view the throat clearly. The lamp is detachable and replaceable.

There are three blades of different sizes that can be detached and attached to the handle according to your professional requirement. The sizes of the blades are 1, 2 and 3.