



Email : sales@naugramedical.com  
Phone: 9896600003

**Product Name :**  
Stethoscope Single Headed

**Product Code :**  
NMDIG10007



**Description :**

Stethoscope Single Headed: Diaphragm only aluminum chest-piece. Silver head. Adjustable chrome-plated binaurals. Flexible PVC Y tubing. White plastic ear-tips. Latex free. Available color of PVC Y-tubing: Red, Green, Gold, Blue, Purple, Black, Pink, Grey, White, etc. Stethoscope Single Headed, Diagnostic Equipments Exporters, Diagnostic Equipments Manufacturers



**Medical Hospital Equipments Manufacturers**

Website: [SITE](#)