

Product Name :
Emergency and Recovery Trolley**Product Code :**
NMEGY10012**Description :**

Emergency and Recovery Trolley Size: 210L x 65W x 65-95H cm. Height adjusted by foot operated hydraulic pump from 65 cm to 95 cm Frame work made of ERW tubular/ rectangular pipes mounted on heavy duty castors of 15 cm Dia (two with brakes) Gas spring assisted Trendelenburg/Reverse Trendelenburg positions. X-Ray Translucent removable Stretcher, backrest raised on ratchet. Complete with corner buffers, synthetic rubber covered handles, sliding X-ray cassette, accessories tray, oxygen cylinder holder. Telescopic saline stand with Stainless Steel drip rod. Swingway Stainless Steel side rails