



Email : sales@naugramedical.com
Phone: 9896600003

Product Name :
Stretcher On Trolley

Product Code :
NMEGY10017



Description :

Stretcher On Trolley Size : 210L x 56W x 82H cms. Frame work made of ERW tubular pipe of 31 mm. dia mounted on four 15 cms. dia. castors (2 with brakes). Removable stretcher top made of S.S. Sheet pressed on tubular frame. Brackets on alternate ends for IV drip rod. Knockdown construction



Medical Hospital Equipments Manufacturers

Website: SITE