



Email : sales@naugramedical.com
Phone: 9896600003

Product Name :
Suction Unit – Eurovac A (Battery Operated)

Product Code :
NMEGY10022



Description :

Suction Unit – Eurovac A (Battery Operated) Housing: MS Powder Coated Capacity: -600 mm Hg \pm 10 , 14 Ltrs/ min Pump Type: Piston Pump Jars: 1 x 0.5 Ltr Polycarbonate Jar Filter: Bacterial Filter Autoclavable/Reusable Tubing: 8 mm ID x 2 mtr. (Silicon) Vacuum Gauge: 5 cm dia 0-760 mm Hg calibration. Power: 220 V AC, 50 / 60 Hz, 55Watt , 12 V DC (Internal Rechargeable Battery 12V 4.5Amp) Battery: Lead Acid (Ni-cd Optional) Noise Level: 55 dB A \pm 3 Dimension & Weight: 37 x 14 x 30 cms, 5.5 Kg. Overflow Safety: Mechanical / Electronic Acoustic Alarm Accessories: Car Charger & Power Adaptor



Medical Hospital Equipments Manufacturers

Website: SITE