



Email : sales@naugramedical.com
Phone: 9896600003

Product Name :
Oxygen Concentrators and Cylinder

Product Code :
Surgical Instruments0047



Description :

Oxygen Concentrators and Cylinder Oxygen concentration: 90% Oxygen flow: 1~5L/min Output pressure: 30-70kPa Overrunning pressure alarm Noise: 52dB(A) Power supply: AC220V, 50Hz Input power: 500VA



Medical Hospital Equipments Manufacturers

Website: Medical Hospital Equipments Manufacturers,