

Product Name :
Oxygen Concentrator with Double Flow

Product Code :
Surgical Instruments0048



Description :

Oxygen Concentrator with Double Flow Oxygen concentration: ~90% Oxygen flow: 1~5L/min Output pressure: 30-70kPa Overrunning pressure alarm Noise: ~60dB(A) Power supply: AC220V±22V, 50Hz±1Hz Input power: ~500VA

NAUGRA[®]

Medical Hospital Equipments Manufacturers

Website:Medical Hospital Equipments Manufacturers,