

**Product Name :**  
Oxygen Concentrators**Product Code :**  
Surgical Instruments0049**Description :**

Oxygen Concentrators It apply zeolite molecular sieve as shift, adopt advanced physical theory(PSA) to separate oxygen and nitrogen in air. At the same time filter the bad substance in air, so as to get pure oxygen for medical use Oxygen Concentration : 93%±3%. Oxygen flux: 1~3L/min. Output pressure: 0.04~0.05Mpa. Overrunning pressure alarm. Operation noise: 49dbA. Power supply: AC220V±20V, 50Hz±1 Hz. Input power: 350VA

**NAUGRA<sup>®</sup>****Medical Hospital Equipments Manufacturers**

Website: Medical Hospital Equipments Manufacturers,