



Email : sales@naugramedical.com  
Phone: 9896600003

**Product Name :**  
Examination Lamp/Light Floor Stand

**Product Code :**  
Surgical Instruments0071



**Description :**

Examination Lamp/Light Floor Stand Approx. : 100 cm. high, floor stand type. Arm : approx. 75 cm. articulated, spring loaded arm. Voltage : 110/220V. Bulb : max. 100w incandescent bulb, Pin Type. Shade : Approx. 15cm dia, light opening. Color : Off-white. Free cord: approx. 2 yard



**Medical Hospital Equipments Manufacturers**

**Website:**Medical Hospital Equipments Manufacturers,