

**Product Name :**  
C-Arm Compatible Fully Electromatic O.T. Table With  
Accentric Base

**Product Code :**  
NMCBS112005



**Description :**

C-Arm Compatible Fully Electromatic O.T. Table With Accentric Base Technical Constraints Top Dimension : L 1930 xW 553mm Height Adjustmen : 762mm- 1012mm Trendelenburg / Reverse :  $30^{\circ}\frac{1}{2}$  /  $25^{\circ}\frac{1}{2}$  Lateral Tilt :  $20^{\circ}\frac{1}{2}$  /  $20^{\circ}\frac{1}{2}$  Kidney Elevator : 150mm Back Rest (up / down) :  $80^{\circ}\frac{1}{2}$  /  $25^{\circ}\frac{1}{2}$  Leg Rest (up / down) :  $15^{\circ}\frac{1}{2}$  /  $90^{\circ}\frac{1}{2}$  Head Rest (up / down) :  $20^{\circ}\frac{1}{2}$  /  $60^{\circ}\frac{1}{2}$  Power Supply : 24V DC Battery Backup (Optional at extra price) : 2 - 3 hours Standard Accessories : Anaesthetic screen : 1 pc. Shoulder supports : 2 pcs. Lateral supports : 2 pcs. Arm rest (S.S.top) : 1 pc. Knee crutches : 2 pcs Overview : Electromatic O.T. Table is equipped with imported electric actuators for smooth & effortless working of the table. It ensures the actuators working in a constant speed without any noise & is designed to run even at low voltage to carry on the surgery Key Structure: The table has been made from acid proof, stainless steel (the face of the table is easy to clean and immune to disinfection agents) Strong & rust free base of stainless steel is provided to make the working conditions pleasant for the doctors & nurses. The various positions such as Up - Down movements, Trendelenburg - Reverse Trendelenburg, Lateral Tilts, Flex - Reflex, Chair Position as realized by remote control