



Email : sales@naugramedical.com
Phone: 9896600003

Product Name :
C-Arm Compatible Fully Electromatic O.T. Table With
Accentric Base

Product Code :
NMCBS112005



Description :

C-Arm Compatible Fully Electromatic O.T. Table With Accentric Base Technical Constraints Top Dimension : L 1930 xW 553mm Height Adjustmen : 762mm- 1012mm Trendelenburg / Reverse : $30^{\circ} \frac{1}{2}$ / $25^{\circ} \frac{1}{2}$ Lateral Tilt : $20^{\circ} \frac{1}{2}$ / $20^{\circ} \frac{1}{2}$ Kidney Elevator : 150mm Back Rest (up / down) : $80^{\circ} \frac{1}{2}$ / $25^{\circ} \frac{1}{2}$ Leg Rest (up / down) : $15^{\circ} \frac{1}{2}$ / $90^{\circ} \frac{1}{2}$ Head Rest (up / down) : $20^{\circ} \frac{1}{2}$ / $60^{\circ} \frac{1}{2}$ Power Supply : 24V DC Battery Backup (Optional at extra price) : 2 - 3 hours Standard Accessories : Anaesthetic screen : 1 pc. Shoulder supports : 2 pcs. Lateral supports : 2 pcs. Arm rest (S.S.top) : 1 pc. Knee crutches : 2 pcs Overview : Electromatic O.T. Table is equipped with imported electric actuators for smooth & effortless working of the table. It ensures the actuators working in a constant speed without any noise & is designed to run even at low voltage to carry on the surgery Key Structure: The table has been made from acid proof, stainless steel (the face of the table is easy to clean and immune to disinfection agents) Strong & rust free base of stainless steel is provided to make the working conditions pleasant for the doctors & nurses. The various positions such as Up - Down movements, Trendelenburg - Reverse Trendelenburg, Lateral Tilts, Flex - Reflex, Chair Position as realized by remote control



Medical Hospital Equipments Manufacturers

Website: Medical Hospital Equipments Manufacturers,