

Product Name :
Electric Orthopedic Imaging Carbon Fiber Operation Table

Product Code :
NMCBS112014



Description :

Electric Orthopedic Imaging Carbon Fiber Operation Table Technical specifications Panel material : Carbon fiber Table Top (Length X Width) : 2020 $\frac{1}{2}$ 10 mmX535 $\frac{1}{2}$ 10mm Height Adjustment : 995 $\frac{1}{2}$ 20mm- 695 $\frac{1}{2}$ 20mm Height range : 300 $\frac{1}{2}$ 20 mm(electric) Reversed Trendelenburg : 30 $\frac{1}{2}$ 1 $\frac{1}{2}$ 2 $\frac{1}{2}$ (electric) Trendelenburg : 30 $\frac{1}{2}$ 1 $\frac{1}{2}$ 2 $\frac{1}{2}$ (electric) Lateral Tilt : 25 $\frac{1}{2}$ 1 $\frac{1}{2}$ 2 $\frac{1}{2}$ (electric) Head Section up/down : 45 $\frac{1}{2}$ (manual)/ 90 $\frac{1}{2}$ (manual) Back Section up/down : 65 $\frac{1}{2}$ 1 $\frac{1}{2}$ 2 $\frac{1}{2}$ (electric)/ 20 $\frac{1}{2}$ 1 $\frac{1}{2}$ 2 (electric) Leg Section up/down : 35 $\frac{1}{2}$ 1 $\frac{1}{2}$ 5 $\frac{1}{2}$ (manual)/ 90 $\frac{1}{2}$ (manual) Translation distance : 330mm (electric) Bulid-in battery : 24V 15AH Kidney Section Elevation : 120mm Supply : AC 220 $\frac{1}{2}$ 10%,50HZ

NAUGRA[®]

Medical Hospital Equipments Manufacturers

Website:Medical Hospital Equipments Manufacturers,