

**Product Name :**  
Electrically Operated OT Table with Double Layer Table  
Top Suitable for X-Ray Examintaion

**Product Code :**  
NMCBS112016



**Description :**

Electrically Operated OT Table with Double Layer Table Top Suitable for X-Ray Examintaion Standard Features: Suitable for surgeries of head, neck, pleuroperitoneal cavity, perineum, limbs, gynecology and obstetrics department, urology thoracic surgery etc. Double layer table top suitable for x-ray examination. Framework & accessories made of stainless steel. Mattress is resistant against acids, alkaline solutions and is anti-bacterial & can be cleaned easily. Various positions like height adjustment, trendelenburg, reverse trendelenburg, lateral tilt, back rest all controlled by linear power device. Equipped with auto chargeable device & battery back up. Leg sections are detachable. In leg section, leg plates are dismountable and separable Technical specifications Table Top (Length x Width) : 2050mm x 480mm Height Adjustable : 700mm-950mm Trendelenburg :  $\pm 25^\circ$  Reverse Trendelenburg :  $\pm 25^\circ$  Lateral Tilt :  $\pm 20^\circ$  /  $\pm 20^\circ$  Back Rest (up/down) :  $\pm 80^\circ$  /  $\pm 15^\circ$  Head Rest (up/down) :  $\pm 80^\circ$  /  $\pm 80^\circ$  Kidney Bridge Elevation :  $\pm 120$ mm Leg Rest adjustment :  $90^\circ$  Power : 220V / 50HZ Foot rest can be kept horizontally or perpendicularly to the Leg Section. Standard Accessories : Knee Crutch : 2 pcs Body Strap : 2 pcs Arm Rest : 2 pcs Anesthesia Rack : 1 pc. Controller : 1 pc