



Email : [sales@naugramedical.com](mailto:sales@naugramedical.com)  
Phone: 9896600003

**Product Name :**  
Sphygmomanometer Aneroid Mobile

**Product Code :**  
NMDIG10016



**Description :**

Sphygmomanometer Aneroid Square Shaped: Mobile, ABS Round dial face adjustable for easy reading, stand on aluminum tube with 5 plastic castors, metal basket, standard latex bulb, inflation valve with spring, standard end valve, PVC tubing. Sphygmomanometer Aneroid Mobile, Diagnostic Equipments Exporters, Diagnostic Equipments Manufacturers



**Medical Hospital Equipments Manufacturers**

Website: [SITE](http://SITE)