



Email : sales@naugramedical.com

Phone: 9896600003

Product Name :
Sphygmomanometer Aneroid Mobile

Product Code :
NMDIG10016



Description :

Sphygmomanometer Aneroid Square Shaped: Mobile, ABS Round dial face adjustable for easy reading, stand on aluminum tube with 5 plastic castors, metal basket, standard latex bulb, inflation valve with spring, standard end valve, PVC tubing. Sphygmomanometer Aneroid Mobile, Diagnostic Equipments Exporters, Diagnostic Equipments Manufacturers



Medical Hospital Equipments Manufacturers

Website:Medical Hospital Equipments Manufacturers,