

**Product Name :**  
Electro Hydraulic Operated Multi Purpose OT Table

**Product Code :**  
NMCBS112022



#### Description :

Electro Hydraulic Operated Multi Purpose OT Table Standard Features: Suitable for thoracic surgery, abdominal surgery, orthopaedics, cerebral surgery, ophthalmology, otorhinolaryngology, obstetrics and gynecology and urinary surgery. X-ray permeable table top cassette guide rail available under the table top. Framework & accessories made of stainless steel. Mattress is resistant against acids, alkaline solutions and is anti-bacterial & can be cleaned easily. Various positions like height adjustment, trendelenburg, reverse trendelenburg, lateral tilt, back rest all controlled by electric hydropneumatic device. Head & Leg sections are detachable. In leg section, separable & dismountable leg plates are optional (at extra cost) Technical specifications Table Top (Length x Width) : 1960mm x 480mm Height Adjustable : 690mm-930mm Trendelenburg :  $\pm 25^{\circ}$  Reverse Trendelenburg :  $\pm 20^{\circ}$  Lateral Tilt :  $\pm 18^{\circ}$  /  $\pm 18^{\circ}$  Back Rest (up/down) :  $\pm 75^{\circ}$  /  $\pm 10^{\circ}$  Head Rest (up/down) :  $\pm 45^{\circ}$  /  $\pm 50^{\circ}$  Kidney Bridge Elevation : 80mm Leg Rest adjustment :  $\pm 90^{\circ}$  Detachable Standard Accessories : Knee Crutch : 2 pcs Body Strap : 2 pcs Arm Rest : 2 pcs Anesthesia Rack : 1 pc. Controller : 1 pcs