

Product Name :
Sphygmomanometer Aneroid Single Tube**Product Code :**
NMDMP10013**Description :**

Sphygmomanometer Aneroid Single Tube: 65mm Dia. dial. Black Powder Coated Large Handle. Std. Bladder Single Tube, Std. Latex or PVC Bulb. Vinyl zipper case, standard end valve. Sphygmomanometer Aneroid Single Tube, Diagnostic Equipments Exporters, Diagnostic Equipments Manufacturers

NAUGRA[®]**Medical Hospital Equipments Manufacturers**Website: SITE