

Product Name :
Range of Motion Knee Brace**Product Code :**
NMBSK1012**Description :**

Range of Motion Knee Brace Size : Universal Feature : Dial screw to give desired angle for range of motion
Foam padded splints to give comfortable fitting Application : Pre and post operative and rehabilitation care
Anterior cruciate ligament (ACL) and Posterior cruciate ligament tear (PCL) Used for Total Knee Replacement (TKR) rehabilitation care

NAUGRA[®]**Medical Hospital Equipments Manufacturers**Website: MedicalHospitalEquipmentsManufacturers.com