

**Product Name :**  
Knee Immobilizer**Product Code :**  
NMBSK1014**Description :**

Knee Immobilizer Size : S, M, L, XL, XXL Feature : High quality foam, fused fabric with metallic splinting for comfort and rigid support Velcro closures for easy and quick application Anatomically shaped splinting for proper fit. Application : An excellent product for support and stability to the injured knee Emergency care in case of fracture, muscle / ligament /tendon injury or dislocation of the knee Post operative rehabilitative care

**NAUGRA<sup>®</sup>****Medical Hospital Equipments Manufacturers**Website: [SITE](http://SITE)