



Email : sales@naugramedical.com
Phone: 9896600003

Product Name :
U SPLINT

Product Code :
NMBSK1018



Description :

U SPLINT Size: Short arm, Long arm, Short leg, Long leg. Features: Malleable, semi-rigid, light metal splint provides support in favourable anatomical position. Provides immobilization to the injured part. Splint is X - ray translucent. Applications: Emergency splinting of dislocated / fractured upper or lower limbs



Medical Hospital Equipments Manufacturers

Website:Medical Hospital Equipments Manufacturers,