

Product Name :
Finger Cot**Product Code :**
NMABW10002**Description :**

Finger Cot Size : S, M, L Feature : V shaped malleable Aluminium splints designed for easy application Poly foam padding for comfort Fits any finger, meant to keep straight for given period of time. Application : Fractures of distal inter phalangeal joints Collateral ligament injury of the distal inter phalangeal joint. Protection from Nail bed injuries of fingers

NAUGRA[®]**Medical Hospital Equipments Manufacturers**Website: MedicalHospitalEquipmentsManufacturers.com