



Email : sales@naugramedical.com
Phone: 9896600003

Product Name :
Frog Splint

Product Code :
NMABW10003



Description :

Frog Splint Size: S, M, L Feature: Foam padded aluminium splint in unique shape to cover the dorsal surface of finger to remain open Four extended arms to ensure good fitting. Application: Hyper flexion injuries Swan neck deformity To hold distal inter phalangeal joints in correct alignment Collateral ligament injuries of distal inter phalangeal joint



Medical Hospital Equipments Manufacturers

Website: SITE