

Product Name :
Finger Extension Splint**Product Code :**
NMABW10005**Description :**

Finger Extension Splint Size : Universal Feature : Malleable aluminium used to keep interphalangeal joints in right position Polyfoam padding for comfort Fits any finger Can be moulded by bare hands Application : Extensor tendon injuries of the finger Collateral ligament injuries of the inter phalangeal joint Boutonniere's deformity Immobilization to the fractured finger

NAUGRA[®]**Medical Hospital Equipments Manufacturers**Website: MedicalHospitalEquipmentsManufacturers.com