



Email : sales@naugramedical.com  
Phone: 9896600003

**Product Name :**  
Macintosh Laryngoscope Blades

**Product Code :**  
NMLYG10B02



**Description :**

Macintosh Type Curved Laryngoscope Blade Dull Finish Size 0, 1, 2, 3, 4, 5 Macintosh Laryngoscope Blades, Macintosh Laryngoscope Blades Exporters, Macintosh Laryngoscope Blades Manufacturers

**Technical Specification :**

Macintosh Laryngoscope Blades combine durable stainless steel with large diameter, high performance full integrated fiber optic bundles.

Designed with sterilization in mind, integrated system minimizes fluid entrapment. Epoxy sealed channel is resistant to temperatures up to 392°F.

Sterilize in Steam, EO and Gas Plasma Systems

Green system compliant

Reinforced ball-lock system

Flawless pinhole free Ni/Cr plating

Select quality 304 Stainless Steel



**Medical Hospital Equipments Manufacturers**

Website: [SITE](http://SITE)