



Email : sales@naugramedical.com  
Phone: 9896600003

**Product Name :**  
Miller Laryngoscope Blades

**Product Code :**  
NMLYG10A02



**Description :**

Miller Type Straight Laryngoscope Blades, Polished Finish - Size 0, 1, 2, 3, 4 Miller Laryngoscope Blades, Macintosh Laryngoscope Blades Exporters, Macintosh Laryngoscope Blades Manufacturers

**Technical Specification :**

The blades of the laryngoscope are made of 304-grade stainless steel. The material makes the blade safe to use and durable.

The handle is made of durable stainless steel. It sports tiny grooves all over its body to allow a good grip.

An LED lamp is fitted along the blade to help a trained medical practitioner view the throat clearly. The lamp is detachable and replaceable.

There are three blades of different sizes that can be detached and attached to the handle according to your professional requirement. The sizes of the blades are 1, 2, 3 and 4.



**Medical Hospital Equipments Manufacturers**

Website: [SITE](#)