

Product Name :
ULTRASOUND THERAPY UNIT**Product Code :**
PEETE10006**Description :**

ULTRASOUND THERAPY UNIT The AutoSound automatically moves the ultrasound output with multiple crystals, Unit with a patented design allows 1MHz, 3MHz, outputs as a Hands-Free therapy. SPECIFICATIONS: Strap the applicator in place for Hands-free Ultrasound at normal therapeutic level. The unique 1 & 3 MHz frequency chain treatment allows you to sweep the entire tissue with deep and superficial therapeutic ultrasound in a single treatment. The Auto Sound automatically moves the Ultrasound output at the ideal 4cm/second (never too slow nor too fast) plus its unidirectional treatment angles. Fitted with patient shut- off control for added safety

NAUGRA[®]**Medical Hospital Equipments Manufacturers**Website: [MedicalHospitalEquipmentsManufacturers,](http://MedicalHospitalEquipmentsManufacturers.com)