

Product Name :
NIRMAL ANKLE AND LEG EXERCISER**Product Code :**
PEEEE10002**Description :**

NIRMAL ANKLE AND LEG EXERCISER Aluminium moulded special designed unit, mounted on cast iron base. Unit is adjustable in height and can permit a full range of motion in eversion, inversion, dorsiflexion and plantar flexion. Resistance can be applied to these motions by addition of weights to the weight pans. Comes with a set of weights (6kg.) in disc form. Aluminium parts are bright polished and other parts finished in stove enamel

NAUGRA[®]**Medical Hospital Equipments Manufacturers****Website:**Medical Hospital Equipments Manufacturers,