



Email : sales@naugramedical.com  
Phone: 9896600003

**Product Name :**  
Laryngoscope Battery Handles Penlight

**Product Code :**  
NMLYG10H02



**Description :**

Battery Handles Penlight Size. Suitable for Laryngoscope Blades. Polished Finish made of Brass with Chrome Plated OR Stainless Steel Dull Finish. Laryngoscope Battery Handles Penlight, Laryngoscope Handles Exporters, Laryngoscope Handles Manufacturers

**Technical Specification :**

Lightweight and convenient Laryngoscope Handle with Xenon Halogen Lamp or LED illumination. Single-piece, stainless steel blade ensures integrity and minimizes corrosion with Knurled handle for a secure grip.



**Medical Hospital Equipments Manufacturers**

Website: Medical Hospital Equipments Manufacturers,