

Product Name :
ULTRA VIOLET LAMP**Product Code :**
PEEHC10005**Description :**

ULTRA VIOLET LAMP Floor stand model for Ultraviolet radiation treatments. UV bulb of 300w is fitted in metal casing with auto-matic cut off timer of 15 minutes. Adjustable arm permits wide range of adjustments. Mounted on heavy duty base fitted with ball bearing casters. Comes with suitable power cord

NAUGRA[®]**Medical Hospital Equipments Manufacturers**Website: SITE