



Email : sales@naugramedical.com
Phone: 9896600003

Product Name :
MOIST HEAT THERAPY UNIT

Product Code :
PEEHC10007



Description :

MOIST HEAT THERAPY UNIT Made out of heavy gauge stainless steel sheet, double walled, fully insulated in between. Overall size 45cm x 30cm x 75cm high. Fitted with 2,000 watts immersion heater, pilot light and thermostat for heat control. Comes with 8 stream packs. (Five standard, one half, one neck and one full back size)



Medical Hospital Equipments Manufacturers

Website:Medical Hospital Equipments Manufacturers,