



Email : sales@naugramedical.com  
Phone: 9896600003

**Product Name :**  
MOIST HEAT THERAPY UNIT

**Product Code :**  
PEEHC10007



**Description :**

MOIST HEAT THERAPY UNIT Made out of heavy gauge stainless steel sheet, double walled, fully insulated in between. Overall size 45cm x 30cm x 75cm high. Fitted with 2,000 watts immersion heater, pilot light and thermostat for heat control. Comes with 8 stream packs. (Five standard, one half, one neck and one full back size



**Medical Hospital Equipments Manufacturers**

Website: SITE,