

Product Name :
MOIST HEAT THERAPY UNIT**Product Code :**
PEEHC10009**Description :**

MOIST HEAT THERAPY UNIT Made out of heavy gauge stainless steel sheet, double walled, well insulated in between. Overall size 58cm x 40cm x 76cm high. Fitted with 2,000 watts immersion heater, pilot light and thermostat for heat control. Comes with 12 stream packs of different sizes (6 standard size, 2 half size, 2 neck size and one full back size). SPECIFICATIONS: Dimensions: Over all, approx.60cm x40cm x80cm. Heat Control: Thermostate to control Temprature from 30 to 110 deg.Cen.with 2 Pilot Lamps. Heater: 2 kw. specifically designed Heater. Hot Packs: Total 12 of different sizes as a) Six, Standard size (25x30cm) b) Two, Half Size (12x30cm) c) Two, Neck Contour (60x15cm) d) Two, Full Back (25x50cm) Accessories: Two forceps to get-out Hot-packs from hot water & Power Pin Plug. Finish: Polished finish inside

NAUGRA[®]**Medical Hospital Equipments Manufacturers**

Website:Medical Hospital Equipments Manufacturers,