

Product Name :
PARAFFIN WAX BATH**Product Code :**
PEEHC10011**Description :**

PARAFFIN WAX BATH Internal wax tank made out of stainless steel sheet. The top has a aluminum cover and the edges of the bath has a laminated wooden rim, works on electricity 220-240 V. Provided with special wax immersion type heater and thermostat to control temperature. Externally finished oven baked. Wax Chamber size 40cm x 20cm x 20cm., Comes with Three kg wax

NAUGRA[®]**Medical Hospital Equipments Manufacturers**Website: [MedicalHospitalEquipmentsManufacturers,](http://MedicalHospitalEquipmentsManufacturers.com)