

Product Name :
AXIAL SHOULDER EXERCISER**Product Code :**
PEEEE10031**Description :**

AXIAL SHOULDER EXERCISER For shoulder & supination-pronation exercises. Used in the same way as a conventional shoulder wheel & wrist machine. A wall mounting unit suitable for small clinics is adjustable in height & radius and fitted with calibrated sensitive resistance mechanism, resistance is controllable from zero to maximum by tightening the strap. A 360 deg. scale is provided on the drum to measure and record the degree of revolution from either direction. Adjustment bars are chrome plated. Comes with wall fixing plates

NAUGRA[®]**Medical Hospital Equipments Manufacturers**Website: MedicalHospitalEquipmentsManufacturers.com