

Product Name :
WRIST CIRCUMDUCTOR**Product Code :**
PEEEE10041**Description :**

WRIST CIRCUMDUCTOR Cast iron wheel of 30cm diameter with drum stand is mounted on a laminated wooden base size:65cm x 22cm. Has a small adjustable handle fitted to the wheel. Knob is provided on the drum to control resistance from zero to maximum. A hollowed foam padded platform size: 22cm x 13cm x 15cm high is fitted with two straps to hold arm and is fitted over the base in front of the wheel for wrist circumduction. Metal parts oven baked finish

NAUGRA[®]**Medical Hospital Equipments Manufacturers**

Website:Medical Hospital Equipments Manufacturers,