

Product Name :
HYDROTHERAPY TANK**Product Code :**
PEEHT10002**Description :**

HYDROTHERAPY TANK Full body hydrotherapy tank is fabricated out of 1.6 mm thick stainless steel sheet, tank is electrically welded and shape is modified figure "8". the overall dimensions are 264cm x 195cm x 86cm high and inside dimension are length to inset 218cm, width upper part 180cm, width of middle narrow part 89cm, width of lower part 135cm and depth is 56cm. It is equipped with 2 turbine ejectors and elevators on ball bearing carriages, one digital thermometer, two water inlets, drain and overflow pipes and standard accessories

NAUGRA[®]**Medical Hospital Equipments Manufacturers**Website: MedicalHospitalEquipmentsManufacturers.com