



Email : sales@naugramedical.com  
Phone: 9896600003

**Product Name :**  
WHIRLPOOL BATH

**Product Code :**  
PEEHT10004



**Description :**

WHIRLPOOL BATH Extremity Whirl-pool is long, shallow & accessible to fit a variety of individual needs. Designed to treat the entire body, the unit can be beneficial for those suffering from Hip & Back conditions as well as upper Back, Neck & Shoulder problems. The heavy gauge stainless steel tank size: 120cm x 60cm x 60cm, is fitted with one motorized turbine ejector and aerator, digital temperature indicator cum controller and heater. Mounted on heavy duty 10cm dia wheels. SPECIFICATIONS: Construction: 1.6mm thick Stainless Steel Sheet Tank is continuously welded with Argon Gas Arc welding & mounted on 10cm dia. 4 heavy duty casters. Turbine: Turbine Ejector-Aerator is powered by a vertically mounted 1/2 hp elc. motor & mounted on a spring loaded turbine elevator for height adjust. Size: 120cm x 60cm x 60cm inside the Tank. Water Heater: 2kw. special Heater, with cover. Temperature Control: Digital Temperature Controller cum Thermometer. Drainage: Drainage system for emptying. Accessories: Anodized Aluminium adjustable Seat with Arm rest. Power Supply: 220V AC



**Medical Hospital Equipments Manufacturers**

Website: [SITE](http://SITE)