



Email : sales@naugramedical.com
Phone: 9896600003

Product Name :
Ultra Sonic Scaler (With 3 Tips)

Product Code :
DSUS0002



Description :

Ultra Sonic Scaler (With 3 Tips) Special self- contained water supply system makes it portable. Enable to use dental solutions for different treatment. Multifunction for general scaling, Endodontic, Periodontal, and restorative treatments. Solution bottle can be change within 3 seconds. Compatible with EMS tips. Built- in water-cooling system extend PCB lifetime. Auto-boost function provides sufficient power to remove heavy calculus. Top design: all surfaces are easy for clean and maintenance. Autoclavable hand piece can be autoclaved up to 135½C. No Compressor & No additional water system required



Medical Hospital Equipments Manufacturers

Website: Medical Hospital Equipments Manufacturers,