



Email : sales@naugramedical.com
Phone: 9896600003

Product Name :
Light Cure Unit With Cordless Digital Lcd Display (LED based)

Product Code :
DLCU0002



Description :

Light Cure Unit With Cordless Digital Lcd Display (LED based) Curing hand piece. Cordless design, LCD display. The special high performance diode mechanism can polymerize an average of 3mm thickness of material within 10 seconds. The smart three setting: Full Power, Soft-up and Impluse Luminosity – 800-1200 mW/cm². LED Lamp: 5 Watt. LED lamp life time up to 100, 000hrs. Wavelength: 450-490nm. Cures in 10 sec./3mm. Additional Lithium-ion Battery



Medical Hospital Equipments Manufacturers

Website: [SITE](#)