

Product Name :
CPAP with Manual Ventilation

Product Code :
ICU0003



Description :

CPAP with Manual Ventilation

Technical Specification :

7" TFT LCD screen displays the ventilation parameters, alarm information and oscillograms.
Intergrated breathing circuit design, ensure easy operating and keep hygiene.
Multiple parameters monitoring interface ensures every parameter is clear, let operator know the patient conditions in all aspects.
Real time pressure-time oscillogram and high precision oxygen concentration detection function included.
Built-in air and oxygen blender, ensure stable oxygen concentration.
Stable and low noise air compressor creates quiet work environment for doctors and patients, able to support 2-3 CPAP machines at the same time.

NAUGRA[®]

Medical Hospital Equipments Manufacturers

Website:Medical Hospital Equipments Manufacturers,