

Product Name :
CPAP without Manual Ventilation**Product Code :**
ICU0004**Description :**

CPAP without Manual Ventilation

Technical Specification :

12 Inch large screen with desirable color schemes
Easy to operate
Support adult and pediatric ventilation mode
Multi-mode ventilation function (VCV, SIMV, PSV , PCV, CPAP , Manual SIGH)
Intelligence operating system
Air Compressor Support (optional)
Self-testing and calibration.

NAUGRA[®]**Medical Hospital Equipments Manufacturers**Website: [SITE](#)