


**Product Name :**  
Multi-Parameter Monitor

**Product Code :**  
ICU0006



#### Description :

Multi-Parameter Monitor High resolution 12.1" color TFT display Lightweight, compact and portable ECG, SPO2, NIBP, RESP, 2-TEMP, PR Optional : 2-IBP, ETCO2, And thermal printer Built-in rechargeable lithium battery ECG waveforms of 7-leads display on the same screen 72-hours graphic and tabular trends of all parameters 72 alarm events of all parameters recall 32 seconds full-disclosure waveform review 500 NIBP measurement data storage and recall Data and waveforms color be adjustable Arrhythmia analysis and S-T segment analysis Suitable for adult paediatric and neonatal patient Standard: ECG, Respiration, Temperature, SPO2, NIBP, Pulse rate Optional: EtCO2, IBP

#### Technical Specification :

10" color TFT LCD screen, High brightness, Wide viewing, 8 waveforms can be displayed on the same screen. lightness and volume are adjustable independently.  
Streamlined outline design with sense of the times & foldable handle.  
Multi-language interface selection: Chinese, English, Spanish, Portuguese, French, Arabic, Russian, Italian, Turkish, etc.  
User-friendly and man-machine dialogue interface with pop-up menus, dialogue boxes, rotary mouse and etc.  
Low-power consumption design with standby mode, rechargeable high-energy built-in battery.  
blood pressure measurement fast and accurately, with venipuncture function.  
Advanced digital technology, accurately measure SpO2 in low perfusion situation.  
Drug calculation, ST-segment detection, pace maker and arrhythmia analysis.  
Display format include Oxy CRG diagram, 7-lead ECG waveform, big-font mode simultaneously ; can store 5 kinds of user-defined display format;  
Sound and light integrated alarm, Alarm parameters can be adjusted to upper or lower.  
Suitable for adults, pediatrics and neonates



## Medical Hospital Equipments Manufacturers

**Website:**Medical Hospital Equipments Manufacturers,