



Email : sales@naugramedical.com  
Phone: 9896600003

**Product Name :**

Ceiling Shadowless OT Light with 2 Domes, One Having 712mm Dia with 7 Reflectors and Another Dome with 550mm Dia Having 4 Reflectors

**Product Code :**

OT Equipment0032

**Description :**

Ceiling Shadowless OT Light with 2 Domes, One Having 712mm Dia with 7 Reflectors and Another Dome with 550mm Dia Having 4 Reflectors Key Structure This Operation reflector Lamp offers the illumination as compared to other standard models. Halogen Lights suitable for Major as well as Minor Surgery, fully spring balanced cool lights. Sterilizable handle has been provided in the middle of the dome to adjust the focus. A Shadowless light is suitable for all sizes of Operation Theatres as it gives maximum light output with minimum heat dissipation. Parabolic reflectors & color correcting filters. Technical specifications Dome head : 712mm / 550mm Halogen Bulbs : 24V x 70W x 11 nos. Light Intensity at 1 mt : 90,000 Lux / 50,000 Lux Height Adjustment : 714 mm / 914 mm Parabolic Reflector : 180 mm Color Temperature : 4200 K 300 Temp rise in field : 3 - 4 C Control Unit : Transformer



**Medical Hospital Equipments Manufacturers**

Website: Medical Hospital Equipments Manufacturers,